**EASYCAP GmbH** Steingrabenstrasse 14 **DE-82211 Herrsching** Germany

Tel +49 (0) 8153 88702-00 Fax +49 (0) 8153 88702-10 www.easycap.de info@easycap.de Delivery Address: Am Anger 5, DE-82237 Woerthsee-Etterschlag

## **Product Information** Forehead Electrode Grid, EasyCap Order Code FEEGRID-22CH-A

## **Features**

- The grid is designed to be placed and fixed on the forehead
- The grid is self adhesive
- The number of electrodes is 24 (11: right side, 11: left side, 2: nose/nasion for reference and ground)
- The electrode material is sintered silver/silverchloride (Ag/AgCl)
- The grid is designed for single use

To connect the grid to an EEG amplifier EasyCap offers dedicated adaptors.



Fig. 1: fEEGrid depicted left shows upper surface, right side shows the surface with electrodes, facing the skin. Self adhesive parts are covered by protection film.

File: fEEGrid-22Ch-A Info EN Status: 2020-11-30 Page: 1 / 1