



**EASYCAP**

*EEG Recording Caps and Related Products*

EASYCAP GmbH  
Am Anger 5  
82237 Wörthsee  
GERMANY

Tel +49 8153 88702-00  
Fax +49 8153 88702-10  
[www.easycap.de](http://www.easycap.de)  
[info@easycap.de](mailto:info@easycap.de)

# **PARTIAL EASYCAP CONSUMABLES PRICELIST**

**effective December 1st, 2021**

**Please note:**

**Prices in this list are  
WITHOUT any VAT, local taxes, fees,  
shipping or processing costs, ...**

**Please trust your local dealer that they will  
always offer you the most favorable price  
possible!**



# EASYCAP

EEG Recording Caps and Related Products

EASYCAP GmbH  
Am Anger 5  
82237 Wörthsee  
GERMANY

Tel +49 8153 88702-00  
Fax +49 8153 88702-10  
[www.easycap.de](http://www.easycap.de)  
[info@easycap.de](mailto:info@easycap.de)

## TABLE OF CONTENTS

### CONSUMABLES & ACCESSORIES

<a href="#">Electrolyte Gels &amp; Pastes</a> .....	3
<a href="#">Syringes &amp; Needles</a> .....	5
<a href="#">Further Consumables</a> .....	6
<a href="#">Input Boards, Connector</a> .....	7
<a href="#">Chest Belts</a> .....	8



## Electrolyte Gels & Pastes

The most effective method to achieve low impedances still is to clean and abrade the skin underneath each electrode, and then to apply an electrolytic gel or paste.

Technically, any conductive gel could be used to bridge the gap between skin and electrode sensor. In practice, different gels are more or less suited for certain EEG applications.

The gels may be

- clear (only conductive) vs. abrasive (conductive AND abrasive)
  - runny vs. viscous
  - salty vs. salt-free
- and have different types of packaging.



<b>Order No.</b>	<b>Description</b>	<b>List Price</b>
V16	SuperVisc 1 kg High-Viscosity Gel for Active Electrodes	59,00 €
V16-10KT	50 pc. SuperVisc (10 gr.) High-Viscosity electrolyte gel in preloaded syringe with safety-Lock, incl. 4 pcs. blunted needles	50,00 €
V16-300	SuperVisc, High-Viscosity Gel for Active Electrodes 300 g cartridge, with adaptor	22,00 €
V17	Abralyt 2000 (1000 gr.) Salt-free, abrasive electrolyte gel	59,00 €
V17-10KT	50 pc. Abralyt 2000 (10 gr.) Salt-free, abrasive electrolyte gel in preloaded syringe with safety-Lock,	50,00 €
V17-300	Abralyt 2000, Salt-free, abrasive electrolyte gel 300 g cartridge, with adaptor	22,00 €
V19	Abralyt HiCl (1000 gr.) High-chloride, abrasive electrolyte gel	59,00 €
V19-10KT	50 pc. Abralyt HiCl (10 gr.) High-chloride, abrasive electrolyte gel in preloaded syringe with safety-Lock, incl. 4 pcs. Blunted Needles	50,00 €
V19-300	Abralyt HiCl, High-chloride, abrasive electrolyte gel 300 g cartridge, with adaptor	22,00 €

.. / .. to be continued next page



# EASYCAP

EEG Recording Caps and Related Products

EASYCAP GmbH  
Am Anger 5  
82237 Wörthsee  
GERMANY

Tel +49 8153 88702-00  
Fax +49 8153 88702-10  
[www.easycap.de](http://www.easycap.de)  
[info@easycap.de](mailto:info@easycap.de)

<b>Order No.</b>	<b>Description</b>	<b>List Price</b>
V14	Lectron II, tube à 250 gr, salt-free electrolyte gel	12,50 €
V15	Lectron III-10, tube à 250 gr, high-chloride electrolyte gel	12,50 €
EDALYT	Isotonic electrolyte gel for GSR/EDA (100 gr)	35,00 €
V11-TUBE	Elefix, conductive paste for EEG testing especially for long term monitoring/sleep studies, slow to dry out. (180 gr.)	15,00 €
V12-N	SAC2, adhesive electrode cream, (100gr)	10,55 €



# EASYCAP

EEG Recording Caps and Related Products

EASYCAP GmbH  
Am Anger 5  
82237 Wörthsee  
GERMANY

Tel +49 8153 88702-00  
Fax +49 8153 88702-10  
[www.easycap.de](http://www.easycap.de)  
[info@easycap.de](mailto:info@easycap.de)

## Syringes & Needles

For filling the gel or paste directly into an electrode, a simple syringe with plain, cone-shaped tip can be used. If a blunted needle is used (e.g. "V39", please see below), a syringe with safety-lock-tip is obligatory ("V35").



<b>Order No.</b>	<b>Description</b>	<b>List Price</b>
V34	Syringe, 10 ml	0,50 €
V34-2	Syringe, 20 ml	0,50 €
V35	Syringe with safety lock, 10 ml	1,00 €
V39	Blunted Needle, metal 19 mm (16G 3/4)	1,80 €
V39-100	Blunted Needle, metal 19 mm (16G 3/4) (box à 100 needles)	96,00 €
V36	Blunted Needle, metal 13mm	1,50 €



# EASYCAP

EEG Recording Caps and Related Products

EASYCAP GmbH  
Am Anger 5  
82237 Wörthsee  
GERMANY

Tel +49 8153 88702-00  
Fax +49 8153 88702-10  
[www.easycap.de](http://www.easycap.de)  
[info@easycap.de](mailto:info@easycap.de)

## Further Consumables

In the following, please find some further consumables.



<b>Order No.</b>	<b>Description</b>	<b>List Price</b>
V52	Washers, inner diameter 4mm, outer diameter 12mm, (100 pieces)*	15,00 €
V52-L	Washers, inner diameter 4mm, outer diameter 19mm, (100 pieces)	15,00 €
V53	Washers, inner diameter 8mm (100 pieces)*	15,00 €
V54	Washers, inner diameter 12mm (100 pieces)*	15,00 €
V40U	Cotton-swabs, non-sterile (100 pieces)	3,00 €
V30	Isopropylalcohol 70% (100 ml)	2,80 €
V94	Perfektan TB (1000 ml), instrument disinfection	39,00 €

\* sizes of washers:





# EASYCAP

EEG Recording Caps and Related Products

EASYCAP GmbH  
 Am Anger 5  
 82237 Wörthsee  
 GERMANY


Tel +49 8153 88702-00  
 Fax +49 8153 88702-10  
[www.easycap.de](http://www.easycap.de)  
[info@easycap.de](mailto:info@easycap.de)

## Input Boards, Connector Adapters

KI input board (38Ch) + KV (38 single sockets, connector to KI input board) serve as quick-connect-interface between single electrodes and most amplifiers (please refer to the picture below, left side).

EIB input boards are used to connect single electrodes to BrainAmp (picture below, right side).

E105 ("squid adapter" like KV) + E80 (flat ribbon cable) can connect BrainCaps to other amplifiers than BrainAmp.



<i>Order No.</i>	<i>Description</i>	<i>List Price</i>
<b>Universal Connectors to any amplifier with single electrode inputs</b>		
KI	38-Channel-Inputbox (for use with KV)	90,00 €
KV	38-Channel-Electrode board adaptor (for use with KI)	120,00 €
<b>BrainAmp Connectors</b>		
EIB32	32-Channel-Electrode Input Box for BrainAmp	180,00 €
EIB64-A	64-Channel-Electrode Input Box for BrainAmp, switchable between single inputs and cap-connector	450,00 €
E105	Universal Converter, 32 Channels, from BrainAmp-connector/EIB (with flat ribbon cable) to 32 + 2 1.5mmtouchproof safety sockets	225,00 €
E80-#	BrainAmp-Flat Ribbon Cable; (# = wildcard for length in cm) prices given for length up to 30cm / length 31-150cm	35,00 / 38,00 €



# EASYCAP

EEG Recording Caps and Related Products

EASYCAP GmbH  
Am Anger 5  
82237 Wörthsee  
GERMANY

Tel +49 8153 88702-00  
Fax +49 8153 88702-10  
[www.easycap.de](http://www.easycap.de)  
[info@easycap.de](mailto:info@easycap.de)

## Chest Belts

Each cap can be fixed by the integrated chin belt alone. As alternative, chest-belts can be used. They give a more comfortable feeling and they hold the cap in place a little bit better.



<b>Order No.</b>	<b>Description</b>	<b>List Price</b>
GS-40-60	Chest belt set XS, extra small, for chest circumference 40 - 60 cm 1 chest belt, 2 chin straps	82,00 €
GS-50-80	Chest belt set S, small, for chest circumference 50 - 80 cm 1 chest belt, 2 chin straps	82,00€
GS-75-105	Chest belt set L, large, for chest circumference 75-105 cm 1 chest belt, 2 chin straps	82,00€
GS-100-130	Chest belt set XL, extra large, for chest circumference 100 - 130 cm 1 chest belt, 2 chin straps	82,00€